

How to Use Prodigy Autocode Diabetes Testing Kit: A Step-by-Step Guide

Managing diabetes requires a regular monitoring of blood glucose levels. The Prodigy Autocode Diabetes Testing Kit, including the Prodigy Auto Code Blood Glucose Meter, is an easy-to-use and reliable tool for tracking blood glucose levels. Here is a detailed step-by-step guide on how to use the Prodigy Autocode Diabetes Testing Kit.



Step 1: Prepare the Testing Kit

The Prodigy Autocode Diabetes Testing Kit includes a blood glucose meter, a lancing device, lancets, and test strips. Before starting the test, ensure the meter is clean and dry, and the test strip is not expired.

Step 2: Wash Your Hands

Use slightly warm water and soap to clean before testing your blood glucose levels. This will help ensure your hands are clean and free of any substances that could affect the test results.

Step 3: Insert a Test Strip

Insert a test strip into the meter. Then, ensure the meter is turned on, and wait until the meter displays the message "Insert Strip."

Step 4: Lance Your Finger

Use the lancing device to prick the side of your finger. Squeeze the finger slightly to obtain a drop of blood.





 **Diabetic
Corner**



 **AutoCode**
Blood Glucose Monitoring System

Step 5: Apply Blood to the Test Strip

Gently touch the end of the test strip to the drop of blood on your finger. The meter will automatically begin to analyze the blood sample.

Step 6: Read the Results

After a few seconds, the meter will display your blood glucose reading. Write the result in a logbook or on your phone or computer to track your blood glucose levels over time.

Step 7: Eject the Test Strip

Then take the test strip, remove it from the meter, and discard it properly.

Step 8: Clean the Meter

Clean the meter with a damp cloth to remove blood or other debris. Then, store the meter and other testing supplies in a cool, dry place.

Using the *Prodigy Autocode Diabetes Testing Kit* is a straightforward process. By following these steps, you can easily monitor your blood glucose levels and manage your diabetes effectively. Consult your healthcare provider for any questions or concerns about managing your diabetes.

It is also very essential to note that the Prodigy Autocode Diabetes Testing Kit is not a substitute that can replace professional medical advice or treatment. While regular monitoring of blood glucose levels is an important part of diabetes management, it is only one aspect of a comprehensive treatment plan. If you come across any queries or concerns about managing your diabetes, consult your healthcare provider.

Conclusion

The Prodigy Autocode Diabetes Testing Kit, including the Prodigy Auto Code Blood Glucose Meter, is a reliable and easy-to-use tool for monitoring blood glucose levels in people with diabetes. By following these steps and tips for using the kit, you can ensure that you get accurate readings and effectively manage your diabetes.



If you are in search of a quality Diabetes Kit then, you can reach out to Diabetic Corner and find Prodigy Autocode Diabetes Testing Kit. You can call us at 951-847-7022 or write to us at sales@daibeticcorner.com.



Contact Us

Diabetic Corner

11801 Pierce Street Turner, Unit 200,

Corona, CA,USA

Zip:- 92505

Email:- sales@diabeticcorner.com

Ph:-19518477022

Website:- <https://www.diabeticcorner.com/>

